



COMPATIBILITY CHART

- ✓ compatible
— not compatible



LLG-Labware Electronic single channel microlitre pipettes IKA-Pen

LLG-Labware Economy 2.0 Pipette Tips	Cat. No.	4723329	4723330	4723331
0.1-10 µL, transparent, pack of 1000	4668775	✓	–	–
0.1-10 µL, transparent, 10 racks of 96	4668776	✓	–	–
0.1-10 µL, XL, transparent, pack of 1000	6491405	✓	–	–
0.1-10 µL, XL, transparent, 10 racks of 96	6491406	✓	–	–
1-200 µL, transparent, pack of 1000	6326227	–	✓	✓
1-200 µL, length 53 mm, yellow, pack of 1000	6524059	–	✓	✓
1-200 µL, length 50 mm, yellow, pack of 1000	4668777	–	✓	✓
1-200 µL, transparent, 10 racks of 96	4668778	–	✓	✓
1-200 µL, yellow, 10 racks of 96	4668779	–	✓	✓
100-1000 µL, length 84 mm, blue, pack of 500	6524060	–	–	–
100-1000 µL, transparent, pack of 1000	4668780	–	–	–
100-1000 µL, transparent, 8 racks of 96	4668782	–	–	–
100-1000 µL, blue, 8 racks of 96	4668783	–	–	–
100-1000 µL, length 71 mm, blue, pack of 1000	4679861	–	–	–
100-5000 µL, transparent, pack of 250	6326975	–	–	–
100-5000 µL, transparent, pack of 250	9409097	–	–	–
1000-10000 µL, transparent, pack of 100	9409098	–	–	–
LLG-Labware Low Retention Pipette Tips				
10 µL, transparent, sterile, 10 racks of 96	4694457	✓	–	–
10 µL, transparent, reload, 10 trays of 96	4694453	✓	–	–
10 µL, XL, transparent, sterile, 10 racks of 96	4694458	✓	–	–
10 µL, XL, transparent, reload, 10 trays of 96	4694454	✓	–	–
200 µL, transparent, sterile, 10 racks of 96	4694459	–	✓	✓
200 µL, transparent, reload, 10 trays of 96	4694455	–	✓	✓
1000 µL, transparent, sterile, 8 racks of 96	4694460	–	–	–
1000 µL, transparent, reload, 5 trays of 96	4694456	–	–	–
LLG-Labware Low Retention 2.0 Filter Tips				
0.1-10 µL, sterile, transparent, 10 racks of 96	4668784	✓	–	–
0.1-10 µL, XL, sterile, transparent, 10 racks of 96	4668785	✓	–	–
1-20 µL, sterile, transparent, 10 racks of 96	4668786	–	✓	✓
1-100 µL, sterile, transparent, 10 racks of 96	4668787	–	✓	✓
1-200 µL, sterile, transparent, 10 racks of 96	4668788	–	✓	✓
100-1000 µL, sterile, transparent, 8 racks of 96	4668789	–	–	–

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been checked.